

Quanteon configuration

Instrument Configuration (Automatically Detected)

Instrument: NovoCyt Quanteon Optical Configuration: 4025 Default [Rename](#) [Create Copy](#)

Note: Please create a copy of current configuration first before customizing the optical configuration.

Par...	Excit...	Detectio...	Default Alias
B525	488	525/45	FITC
B586	488	586/20	EYFP
B615	488	615/20	PI
B667	488	667/30	PerCP
B695	488	695/40	PerCP-Cy5.5
B725	488	725/40	PerCP-eFluor 7...
B780	488	780/60	PE-Cy7 (B)
R667	637	667/30	APC
R695	637	695/40	Alexa Fluor 680
R725	637	725/40	Alexa Fluor 700
R780	637	780/60	APC-Cy7
V445	405	445/45	Pacific Blue
V525	405	525/45	AmCyan
V586	405	586/20	Pacific Orange
V615	405	615/20	Qdot 605
V667	405	667/30	Qdot 655
V695	405	695/40	Qdot 705
V725	405	725/40	BV711
V780	405	780/60	Qdot 800
Y586	561	586/20	PE
Y615	561	615/20	PE-Texas Red
Y667	561	667/30	PE-Cy5
Y695	561	695/40	PE-Cy5.5
Y725	561	725/40	PE-Alexa Fluor...
Y780	561	780/60	PE-Cy7 (Y)

☒ NovoSampler Q installed

[Change Optical Configuration](#) [OK](#) [Cancel](#)